

DEVICE CONNECTION CABLE TYPE: I3

PHYSICAL CONNECTIONS

Ottobock Electric Wrist Rotator (10S17)

CONNECTION:
Coapt cable connection marked 'Close/Digital' must plug into this side '2' of wrist

CONNECTION:
Coapt cable connection marked 'Rotate' must plug into wrist port marked '3'

CONNECTION:
Coapt cable connection marked 'Open/Digital' must plug into this side '2' of wrist

CONNECTION:
Power/battery cable (not supplied by Coapt) must plug into wrist port marked '1'

Motion Control Standard Wrist Rotator (5010045, 5010054, 5010055)

CONNECTION:
Coapt cable connection marked 'Rotate' must plug into wrist port marked 'WRIST'

CONNECTION:
Coapt cable connection marked 'Open/Digital' must plug into wrist port marked 'CH-A'

CONNECTION:
Power/battery cable (not supplied by Coapt) must plug into wrist port marked 'HAND'

CONNECTION:
Coapt cable connection marked 'Close/Digital' must plug into wrist port marked 'CH-B'

SOFTWARE SETTINGS

WRIST SETTINGS:

No manual settings necessary for the Ottobock Standard Wrist Rotator or the Motion Control Standard Wrist Rotator

I-LIMB SETTINGS:

For general i-Limb Digital operation via Coapt, most settings should be set to defaults in the i-Limb app

For proper operation of Gesture Control via Coapt, ensure 'CC' is selected as its trigger in the i-Limb app

